




CONTECH[®]
Instruments Ltd.
AN ISO 9001-2000 COMPANY
Where innovation never stops...

CA SERIES PRECISION BALANCES

0.01g to 5000g



HIGH PRECISION BALANCES

- ISO GLP compliance.
- Standard bi-directional RS-232 interface.
- Multiple weighing units.
- Polyfunction : % weighing, counting, fill mode, GSM etc.

Model	Capacity	Readability	Pan Size
CA-1502	1500g	0.01g	Ø150 mm
CA-2202	2200g	0.01g	
CA-3102	3100g	0.01g	
CA-4102	4100g	0.01g	
CA-5002	5000g	0.01g	